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
fig. D

Contents

Legends	08
Buttons and LEDs	08
Accessories	08
Curing light components	08
Installation	08
Producers prior to equipment use/reutilization	08
Operation	08
Cleaning	08
Disinfection and sterilization.....	09
Maintenance	09
Recommendations, care and warnings.....	09
Procedure in case equipment functioning is altered.....	09
Sensitivity to predictable environmental conditions in normal working situations.....	09
Procedures in case of equipment disposal.....	09
Transportation and storage conditions.....	09
After installation	10
Preparation prior to delivery	10
Dabi Atlante Authorized Service Network.....	10

Symbols	10
Final recommendations	10
Equipment warranty	10
Technical features	11
Troubleshooting.....	11

Presentation

Dear Client,

Welcome!

We are really proud of your preference for our products. Your familiarization with **DB 686** will help maximize its performance and increase its durability and benefits. Therefore, before starting your equipment, be sure to read this manual very attentively.

Important

All the information, illustrations and specifications contained in this manual are based on data existing at the time of its publication. We reserve the right to make any modifications deemed necessary both on the product and in this manual without previous notice.

Legends

Buttons and LEDs (fig. B)

- ① Power button
- ② Green light – power indicator
- ③ Pulse mode: LED blinking steadily
- ④ Gradual mode: light intensity increases gradually
- ⑤ Intense mode: maximum light intensity
- ⑥ Mode selection button

Accessories (fig. C)

- ① Eye protector
- ② Charger

Curing light components (fig. D)

- ① Light guide tip
- ② Handpiece
- ③ Battery
- ④ Charger base
- ⑤ Charge indicator light

Installation

Product installation (fig. D)

Connect the charger to the charger base and to an appropriate power source. Place the charger base on a firm surface.

Screw the battery onto the lower part of the handpiece and insert the light guide tip into the higher part. When the light guide tip is in the correct position, a short clicking sound will be heard.

Put the eye protector into the light guide tip and insert the equipment into the charger base.

- ⚠ Charge the battery during 12 hours before using the equipment for the first time.

Procedures prior to equipment use/reutilization

Prior to equipment use/reutilization, follow the operation and cleaning instructions in this manual.

Operation

Select the desired mode by pressing the mode selection button (fig. B, item 6). The standard mode is the intense mode.

After mode selection, take the handpiece to the patient's mouth and position the light guide tip at a distance of 1-2mm from the material to be cured. Press the power button (Fig. B, item 1) to start operation.

In order to interrupt power supply during the procedure, press the power button again (fig. B, item 1).

After an operation period, the power button can be pressed in order to initiate the next period; however, the equipment must not be used more than 10 sequential times. After such period, wait a few seconds so that the equipment will not overheat.

- ⚠ Keep the equipment out of the charger base until it indicates the necessity of recharge. In this case, the equipment will beep twice every 4 seconds. Even beeping, the equipment will continue to operate but its power will not be ideal for curing.

In case of a long period without use, remove the light guide tip and separate the handpiece and the battery, putting them away under adequate conditions.

Operation Time:

Mode	Time	
	DB 685	DB 686
Pulse mode	20s	20s
Gradual mode	15s	20s
Intense mode	10s	20s

Charger Base

When the equipment is inserted on the charger base, the charge indicator light (fig. D, item 5) will turn on.

If the light is red, it means that the battery is being charged. In case it is green, it indicates that charging has been completed.

Automatic shut off

The equipment will be turned off automatically in case operation is interrupted for 3 minutes. To turn it on again, press the power button.

Battery warning

When the green light power indicator (fig. B, item 2) blinks and 2 beeps are emitted, it means that battery is low. In case that happens during the curing procedure, extend the curing period and recharge the battery in time.

Cleaning

Clean the equipment surface using a clean soft cloth moistened in Aplic Odonto by Dabi Atlante or another product with similar chemical properties.

Do not place the light conduction tip in contact with any rigid item, and after its use check for any material remaining on the terminal. In case any material exists, clean it but do not use metal instruments.

If the light conduction tip is damaged, it must be replaced.

There should be no dirt on the golden surface

of the printed circuit board on both sides of the battery or on the connection chip. Cleaning must be performed by using alcohol.

Disinfection and Sterilization

When decontaminating the handpiece and the battery, they should be completely separated.

The light conduction tip must be cleaned and sterilized at 134°C before it is used on the next patient.

Maintenance

Preventive maintenance

In order to reduce equipment faults and increase durability, call Dabi Atlante Authorized Technical Assistance and set a regular preventive maintenance plan.

-  To guarantee total light intensity set technical visits every 6 months.

Corrective maintenance

In case the equipment presents any faults that are not listed in "Troubleshooting" in this manual, contact Dabi Atlante Authorized Technical Assistance.

-  Do not open the equipment and/or try to repair it as such action may aggravate the problem or even generate further malfunctioning.

Circuit diagrams, list of parts, components, others

In case they are necessary, ask Dabi Atlante for electric diagrams and/or lists of parts, components or other information, which will be supplied by agreement.

Recommendations, care and warnings

This equipment must be used only by properly trained personnel under supervision by qualified dental staff.

It is advised that the dental professional instruct his patient or other individuals who have access to the equipment as regards the caution to be taken during its handling.

This equipment is not appropriate for use in areas at risk of explosion.

Prior to the first use and/or after long work interruptions, such as vacation periods, clean and disinfect the equipment.

Do not directly expose the equipment to sunlight as that will cause early wear.

Prior to connecting the charger base to the power supply, make sure that the main voltage is in accordance with the charger's voltage range.

Do not look at the light without the eye protector as it may be harmful.

Prior to use, wait until the unit reaches room temperature, especially when it is taken from cold conditions into warmer conditions.

Never use a battery other than the original one; otherwise your equipment may be damaged.

Do not place the unit near flammable items, and do not place the light in contact with flammable solutions.

In order to prevent electric shock, do not insert any other parts into the unit.

Do not use the equipment on patients carrying cardiac pacemakers.

Procedure in case equipment functioning is altered

In case the equipment presents any faults, firstly check whether they are listed under "Troubleshooting". If the problem cannot be solved, turn off the equipment, unplug it from the power source outlet and call for Dabi Atlante Authorized Technical Assistance.

Sensitivity to predictable environmental conditions in normal working situations

This equipment is not sensitive to electric, electrostatic or pressure interferences as long as the cleaning, maintenance and transportation instructions in this manual are observed. However, electromagnetic environments may interfere with its normal operation.

Procedures in case of equipment disposal

When the equipment is put out of use, it must be disposed of at an appropriate site (according to local regulations), thus preventing its inadequate use and environmental contamination.

Transportation and storage conditions

The equipment must be transported and stored:

- Carefully in order to prevent falls and impact.
- By keeping the side indicated by the arrow facing upwards.
- By observing maximum piling limits indicated on the packaging.
- Protected against moisture, rain, water sprinkling and wet ground.

Transportation/storage environmental conditions

Transportation and storage environmental temperature range	0°C a + 55°C
Transportation and storage relative humidity range	0% to 90% (non-condensing)
Atmospheric pressure range	500hPa to 1060 hPa (375 mmHg to 795 mmHg)

After installation

Keep the equipment protected from sunlight and rain.

Environmental operation conditions

Environmental temperature functioning range	+15°C to +30°C
Environmental temperature range recommended by Dabi Atlante	+21°C to + 26°C
Functioning relative humidity range	10% to 80% (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

Environmental accommodation operation conditions (between operations)

Environmental temperature accommodation range	+5°C to +45°C
Environmental temperature range recommended by Dabi Atlante	+15°C to +30°C
Functioning relative humidity range	30% to 75% (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

Preparation prior to delivery





This equipment has been revised and approved according to the Dabi Atlante Quality Inspection and Control Program for better performance.

Dabi Atlante Authorized Service Network

Installation and all maintenance service for Dabi Atlante equipment must be given by Dabi Atlante Authorized Technical Assistance, otherwise warranty will be compromised. Use the attached insert or access the website www.dabi.com.br to find the nearest Dabi Atlante Authorized Technical Assistant.

Symbols

Use the icons below to identify the symbols on your equipment:

	"Fragile" Located on the packaging side. It determines that equipment must be carefully transported in order to prevent falls or shocks
	"Protect against moisture" Located on the packaging side. It determines that the equipment must be protected against any type of moisture during transportation and storage
	"This side up" Located on the packaging side. It determines that the packaging must be handled by keeping the side shown by the arrow facing upwards
	"Maximum piling" Located on the packaging side. It determines the maximum number of boxes that can be piled up during transportation and storage
	"Temperature limit" It determines the temperature limits at which the equipment must be stored or transported
	"Maximum Moisture" Determinate the maximum moisture level that the package should be submitted when stored ou transported.

Final recommendations

In order to ensure greater durability for your equipment, use only original Dabi Atlante replacement parts, which are in accordance with the technical standards and specifications required by the manufacturer. Also, use only the services provided by the Dabi Atlante Authorized Technical Assistance Network as it counts on trained assistants and specific tools for the appropriate maintenance of your equipment.

Equipment Warranty

The "Warranty Terms" are found on the back of the certificate accompanying the product. Any alterations or concessions involving the warranty as well as pertaining conditions and authorizations, either oral or written, are definitely not permitted if not previously documented and approved of by Dabi Atlante.

Technical features

	DB 686	DB 685
Classification	Class II (medium risk)	Class II (medium risk)
Charger's voltage range	100-240Vac	100-240Vac
Frequency	50/60Hz	50/60Hz
Power	1.5VA	1.5VA
Lithium battery	3.7V 2000mAh	7.4V 1800mAh
Charging time	5 hours	5 hours
Charge duration	400 cures of 20s*	400 cures of 20s*
Light power source		
Maximum power consumption	3W	5W
Light wavelength	440-480nm	440-480nm
Light intensity	>500mW/cm ^{2**}	>1000mW/cm ^{2**}
Average durability	10.000 hours (if each operation lasts 20s, the equipment can be used up to 1,800,000 times)	10.000 hours (if each operation lasts 20s, the equipment can be used up to 1,800,000 times)
Light guide tip	Ø 08 mm angle: 60°	Ø 08 mm angle: 60°

*Depending of the battery life, it may change

**In the operation conditions indicated in this Manual

Troubleshooting

Faults	Causes	Solutions
Equipment does not operate properly	Low battery	Recharge battery.
	Inadequate contact with battery	Clean the circuit board, the battery connector and the handpiece connector. Connect the handpiece and the battery again.
Charging fails	Inadequate power source contact	Connect the base to the power supply adequately.
	Inadequate contact with battery	Clean the corresponding circuit board and the base connector.

Technical responsibility: Caetano Barros Biagi - CREA-SP: 5061859382D

www.dabi.com.br
Tel.: +55 (16) 3512 1212